



Shenzhen Well Dental Laboratory

Address: 518 Wodu Zonghe Building,Henggang Street,Longgang District,Shenzhen,Guangdong, China,518115

NO.: _____

Dr.'s Name: _____

Patient: _____

American teeth number:

| | | | | | | | | | | | | | | | |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 |
| 32 | 31 | 30 | 29 | 28 | 27 | 26 | 25 | 24 | 23 | 22 | 21 | 20 | 19 | 18 | 17 |

European teeth number:

| | | | | | | | | | | | | | | | |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| 18 | 17 | 16 | 15 | 14 | 13 | 12 | 11 | 21 | 22 | 23 | 24 | 25 | 26 | 27 | 28 |
| 48 | 47 | 46 | 45 | 44 | 43 | 42 | 41 | 31 | 32 | 33 | 34 | 35 | 36 | 37 | 38 |

No lingual metal margin

Shade: _____

If no occlusal clearance:

- reduction preparation
- reduce opposing
- make metal occlusal
- lab should contact us

Fixed restorations:

Porcelain Fused To Metal(PFM): Non-precious(Co-Cr) Semi-precious High precious

Full Metal Cast(FMC): Non-precious(Co-Cr) Semi-precious High precious

All-ceramic: Full milled zirconia Porcelain fused to zirconia EMax

Implant: Cement retained screw retained

Post core

Post with crown/Richmond

Separate post+crown

veneer

onlay

inlay

Temporary crown

maryland bridge

Removable restorations: Upper Lower

Acrylic denture: set up teeth finish

Flexible denture: set up teeth finish

Metal framework: set up teeth finish

Others: Bleaching tray night guard orthodontic appliance

Additional notes:

